

Clinical Concordance for Exclusion of Acute Heart Failure Using the RAMP and Dade Dimension RxL NT-proBNP Assays

M Offner¹, E Winckel¹, A Jung¹, T Lamy², G Arthur³, ML Zucker³

¹Laboratoire de Biochimie Générale et Spécialisée, Hôpital Civil, Strasbourg, France, ²All.DIAG.SA, Strasbourg, France, ³Response Biomedical Corporation, Burnaby, BC, Canada

Introduction

- Accurate evaluation and management of acute dyspnea remains a significant clinical challenge^{1,2}.
- Blood levels of NT-proBNP:
 - are indicative of the degree of heart failure (HF).
 - augment diagnosis and clinical judgment^{3,4}.
 - have been used for risk stratification of patients with acute coronary syndrome and congestive HF^{5,6}.

NT-proBNP (pg/mL) Clinical Cut-offs ^(3,4)			
Age Dependent	<50	50-75	>75
Limit for Probable Acute Heart Failure	>450	>900	>1800
Age Independent Rule Out	<300		

- Compared with central lab testing, point of care testing (POC) allows rapid clinical decision making by reducing the time taken to receive test results.

Aims


- To evaluate clinical concordance of a point of care NT-proBNP assay* (RAMP) with a laboratory assay for NT-proBNP (Dade Dimension RxL).

Methods and Procedures

- Patients (n=91) suspected of HF were enrolled (no age data available).
- EDTA and heparin anticoagulated blood samples were obtained.
- NT-proBNP determined on:
 - Dade Dimension RxL analyzer
 - Heparinized plasma in laboratory
 - RAMP Reader
 - EDTA whole blood (POC)

Results: Method Comparison

- 86 of 91 samples had NT-proBNP levels within the reportable ranges of both test systems (5 - 35,000 pg/mL).



- Linear regression analysis showed excellent correlation (R=0.971) between the two systems with a slope of 1.04.

Results: Clinical Concordance

Clinical concordance analysis of the RAMP NT-proBNP and Dade Dimension RxL results using an age independent rule out value of 300 pg/mL³.

Table 1		Dade Dimension RxL	
		>300 pg/mL	<300 pg/mL
RAMP	>300 pg/mL	70	5
	<300 pg/mL	0	16
		Sensitivity	100.0%
		Specificity	76.2%
		PPV	93.3%
		NPV	100%
		Concordance	94.5%

Clinical concordance analysis showed excellent clinical agreement between the RAMP and the Dimension.

Conclusions

- The results obtained on the RAMP POC system were comparable to those obtained on the Dade Dimension RxL laboratory analyzer.
- The RAMP NT-proBNP whole blood assay is an accurate clinical indicator that can be used to rule-out acute HF in this population.
- With results available in 15 minutes, RAMP NT-proBNP allows for rapid initiation of appropriate patient treatment.

Bibliography

- Mueller C, Scholer A, Laule-Kilian K, et al. Use of B-Type natriuretic peptide in the evaluation and management of acute dyspnea. *N Eng J Med.* 2004;350:647-654.
- Zaphiroy A, Robb S, Murray T, et al. The diagnostic accuracy of plasma BNP and NTproBNP in patients referred from primary care with suspected heart failure: Results of the UK natriuretic peptide study. *Eur J Heart Failure.* 2005;7:537-541.
- Januzzi JL, Camargo CA, Anwaruddin S, et al. The N-terminal Pro-BNP investigation of dyspnea in the emergency department (PRIDE) study. *Am J Cardiol.* 2005;95:948-954.
- Baggish AL, Cameron R, Anwaruddin S, et al. A clinical and biochemical critical pathway for the evaluation of patients with suspected acute congestive heart failure. The proBNP investigation of dyspnea in the emergency department (PRIDE) algorithm. *Crit Path Cardiol.* 2004;3:171-176.
- Jernberg T, Stridsberg M, Venge P, et al. N-terminal pro brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. *J Am Coll Cardiol.* 2002;40:437-45.
- Fisher C, Berry C, Blue L, et al. N-terminal pro B type natriuretic peptide, but not the new putative cardiac hormone relaxin, predicts prognosis in patients with chronic heart failure. *Heart.* 2003;89:879-81.

