



# Evaluation of RAMP Whole Blood Analyzer for Point of Care Troponin I and NT-proBNP Testing



ML Zucker<sup>1</sup>, G Soffiati<sup>2</sup>, G Secco, G Arthur<sup>1</sup>, A Fortunato<sup>2</sup>

<sup>1</sup>Response Biomedical Corporation, Burnaby, BC, Canada. <sup>2</sup>Laboratorio di Chimica Clinica ed Ematologia, Ospedale S.Bortolo, Vicenza, Italy.

## Introduction

- Accurate and rapid measurement of bio-markers aid in the diagnosis of cardiovascular disease<sup>1</sup>.
- This study evaluates the performance of the RAMP Point of Care whole blood tests for Troponin I (TnI), and N-terminal pro-brain natriuretic peptide (NT-proBNP).
- TnI testing methods are not harmonized and normal reference ranges and method comparisons are required to establish clinical cut-offs.
- NT-proBNP assays are harmonized and common clinical cut-offs and ranges have been established<sup>2,3</sup>.

## Methods and Materials

- EDTA and heparinized blood samples were collected from patients for whom testing was clinically indicated, n=46 (TnI) and n=42 (NT-proBNP).
- EDTA samples were collected from subjects without cardiac symptoms to determine the 99<sup>th</sup> percentile cut-off for the RAMP TnI test, n=100.
- Method comparisons and clinical concordance analysis were carried out for RAMP versus:
  - Centaur Ultra Sensitive TnI (Siemens Medical Diagnostic Solutions, Germany).
  - Roche Elecsys 2010 proBNP (Roche Diagnostics, Indianapolis, IN).
- EDTA whole blood samples analyzed on RAMP; heparinized plasma samples analyzed on Centaur (TnI) or Elecsys (NT-proBNP).

## Results: Clinical Agreement


**TnI:** Normal range (99<sup>th</sup> percentile) cut-off on RAMP was determined as 0.16 ng/mL (CV = 18%). A 0.06 ng/mL cut-off was used for the Centaur.

**NT-proBNP:** 300 pg/mL age independent rule out was used on both systems<sup>2</sup>.

	TnI	NT-proBNP
Sensitivity	87.2%	97.1%
Specificity	85.7%	100%
Concordance	87.0%	97.6%

## Results: Troponin I


Passing & Bablok analysis of RAMP versus Centaur for TnI



RAMP TnI correlated strongly with the Centaur TnI (Pearson R=0.981, n=45, one outlier excluded). Strong correlation was maintained when the analysis was confined to samples with TnI levels <5.0 ng/mL by Centaur (Pearson R=0.945; n=31).

## Results: NT-proBNP

Passing & Bablok analysis of RAMP versus Roche Elecsys



RAMP NT-proBNP correlated strongly with the Elecsys proBNP (Pearson R=0.994, n=41, one out of range value excluded) across a wide range of NT-proBNP values. Strong correlation was maintained when analysis was confined to samples <5000 pg/mL (Pearson R=0.987, n=31).

## Conclusions

- The RAMP TnI assay correlates strongly with the Centaur TnI assay and using an appropriately defined cut-off point, provides accurate diagnostic results in less than 20 minutes.
- The RAMP NT-proBNP assay correlates strongly with the Roche Elecsys proBNP assay and provides accurate diagnostic results in 15 minutes.

## Bibliography

- 1 ACEP Clinical Policies, Ann Emerg Med. 2006;48:270-301.
- 2 Baggish AL, Cameron R, Anwaruddin S, et al. A clinical and biochemical critical pathway for the evaluation of patients with suspected acute congestive heart failure. The proBNP investigation of dyspnea in the emergency department (PRIDE) algorithm. Crit Path Cardiol. 2004;3:171-176.
- 3 Januzzi JL, Camargo CA, Anwaruddin S, et al. The N-terminal Pro-BNP investigation of dyspnea in the emergency department (PRIDE) study. Am J Cardiol. 2005;95:948-954.



\*RAMP NT-proBNP Assay: CE marked; Currently under FDA / Health Canada review.